INSTRUCTIONS FOR USE - ENGLISH

IMPORTANT INFORMATION

The words WARNING, CAUTION, and NOTE have special meaning and should be reviewed.

WARNING	A WARNING highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
CAUTION	A CAUTION highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
NOTE	A NOTE supplements and/or clarifies procedural information.

COMPLIANCE STATEMENTS

Medical Device acc. Directive 93/42/EEC

This device is subject to the conformity assessment procedure according to annex II, performed as applicable by the following notified body: DQS Medizinprodukte GmbH (Code 0297)

Personal Protective Equipment acc. Regulation (EU) 2016/425

This PPE has been EC-type examined by the following notified body: CERTOTTICA S.C.R.L., Z.I. Villanova 32013 Longarone (BL) Italy, (Code 0530)

To download the declaration of conformity please visit: www.thigmbh.at

Marking on the sterile barrier system: **THI EN166 3 F**

Marking on the ocular (lens): **THI 1 F**

CODE	DEFINITION
THI	Code of Manufacturer - THI GmbH
EN166	European Norm - Personal Eye Protection Specification - EN 166:2001
3	Field of use - against splashes of liquids
F	Mechanical Strength - resistance against high-speed particles - low energy 45 m/s
1	Optical Class

INDICATIONS FOR USE

The ViVi® Hood HFD is a component of the ViVi® System and is intended to be worn by surgical personnel to provide a barrier between the operating environment and the surgical personnel in order to help protect against contamination and/or exposure of infectious body fluids and harmful microorganisms.

CONTRAINDICATIONS

None known.

LIMITATIONS OF USE

WARNINGS:

- This personal protection system is not intended to be used as a respiratory protection system.
- The lenses of this system provide limited protection to the face from flying particles or the splash of hazardous liquids. If these lenses are worn over ophthalmic spectacles, impacts may be transmitted. The lenses are not shatterproof.
- The hood is not intended for use against oily particulate such as paint mist, oil mist or detergents.
- The hood fabric is not intended to be exposed to open flame and may burn.
- The hood is not intended for use in atmospheres immediately dangerous to life or health.
- The materials contained within this system are not known allergens. However, if an allergic reaction occurs, contact a physician.

FOR USE WITH

The ViVi® Hood HFD is a component of the ViVi® System and is intended to be worn over any ViVi® Helmet.

USER/PATIENT SAFETY

WARNINGS:

- Only trained and experienced healthcare professionals should use this equipment. Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use.
- The risk severity depends on several factors; the end-user shall therefore always evaluate if the PPE meet protective requirements required for the intended activity.
- Scratches and/or other damages to the device may seriously reduce the impact protection level.

- If any part of the device is deemed unfit for use, the whole product should be discarded.
- This PPE shall not be used at extremes of temperatures.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. THI, as a manufacturer, does not recommend surgical procedure or technique.
- Sterile only if package is unopened and undamaged. If contamination is suspected, discard the contaminated hood and use a new hood.
- DO NOT reuse, reprocess, or repack this device. This device is intended for a single use only.
- DO NOT clean or re-sterilize this device. Cleaning or re-sterilization may significantly affect the filtration efficiency and result in tears in the garment or other loss of mechanical integrity. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- Use only THI-approved components and accessories, unless otherwise specified. DO NOT modify any component or accessory.

NOTE:

- If you need more information, contact your local THI Distributor.

SYMBOL DEFINITIONS

The symbols located on the equipment and/or labeling are defined in this section.

SYMBOL	DEFINITION
	Consult instructions for use
	Do not use if package is damaged
	Manufacturer
	Do not reuse
STERILE EO	Sterilized using ethylene oxide
R_x ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
LOT	Batch code
	Use by date
	Hook-and-loop attachment points
QTY	Quantity
REF	Catalog number

INSTRUCTIONS

WARNINGS:

- ALWAYS wear the helmet and hood as a total system to obtain the required personal protection.
- Before and during use, ALWAYS inspect the equipment for rips, tears, cuts, or holes. ALWAYS make sure seams and closing systems are secure and function properly. DO NOT use the equipment if damage is apparent.
- Per hospital protocol, ALWAYS use proper sterile technique.
- Do not contaminate the exterior surface of the hood during donning.
- ALWAYS follow current hospital protocol governing the safe removal, handling, and disposal of potential biohazard waste.

TO DON THE HOOD (FIGURES 1 TO 11)

NOTES:

- The hood can be donned by the user himself. For unpacking the help of an assistant is required.
 - The hood is folded and packaged to permit handling without causing contamination to the exterior surfaces of the hood
- Open the packaging and remove the hood with the help of an assistant. (fig. 1).
 - Point the opening towards you and unfold the wrapping paper. (fig. 2)
 - Roll the sides of the hood to the outside and grab the hood with both hands. (fig. 3)
 - Place the hook opening of the face shield on the hook of the helmet. (fig. 4)
 - Carefully connect the hook-and-loop pieces of the hood onto their corresponding pieces on the helmet. Ensure that you have created a tight fit between the face shield and the helmet. (fig. 5)
 - Carefully pull the front flap of the hood down. (fig. 6)
 - Grab the chest straps and carefully pull the hood with the help of the chest straps over the helmet. (fig. 7)

8. Pull the hood down your back using the chest straps (fig. 8)
9. Straighten up the front flap of the hood (fig. 9)

NOTE:

- Make sure the white filter material of the hood is positioned correctly over the fan intake of the helmet.
10. Tie the chest straps at your chest. Take care that the chest straps pass under your armpits (fig. 10)
 11. Complete the process by donning a surgical gown keeping the hood under the surgical gown. (fig. 11)

STORAGE AND HANDLING

Store the equipment within the specified environmental condition value(s) throughout its useful life. See the Specifications section. To ensure the longevity, performance and safety of this equipment, use the original packaging materials when storing or transporting this equipment.

SPECIFICATIONS

MATERIAL STANDARDS COMPLIANCE	
Barrier Performance (front panel)	• ASTM F1671 (viral penetration) • AATCC 42 (impact penetration)
Flammability	Flammability of Clothing Textiles under United States Consumer Product Safety Commission (US CPSC) 16 Code of Federal Regulations (CFR) Part 1610.
ENVIRONMENTAL CONDITIONS	
Operation and Transportation Temperature Limit	 49 °C
Storage Temperature Limit	 30 °C

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Total Healthcare Innovation GmbH
THI - Total Healthcare Innovation GmbH