CE

R_x ONLY

SURGICAL HELMET

INSTRUCTIONS FOR USE

REF	80100	/ Vivi ® Helmet
REF	80120	/ Vivi ® Helmet w/Headlight
REF	80130	/ Vivi ® Helmet HPL

ENGLISH (EN)
DEUTSCH (DE)
ITALIANO (IT)

Total Healthcare Innovation GmbH

© THI GmbH - 2020/09
REF 20602 | Rev-04

INTRODUCTION

This *Instructions For Use manual* is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians.

Keep and consult this reference manual during the life of the product.

THE FOLLOWING CONVENTIONS ARE USED IN THIS MANUAL:

WARNING	A WARNING highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
CAUTION	A CAUTION highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
NOTE	A NOTE supplements and/or clarifies procedural information.

If additional information, especially safety information, or in-service training is required, contact your distributor.

Trademarks not the property of THI GmbH are the property of their respective owners.

COMPLIANCE STATEMENTS

MEDICAL DEVICE ACC. DIRECTIVE 93/42/EEC

PERSONAL PROTECTIVE EQUIPMENT ACC. REGULATION (EU) 2016/425

In combination with **ViVi® Hood** (REF 80200, 80210, 80300, 80301) or **ViVi® Toga** (REF 80220, 80221, 80222), this device is part of a PPE system. This device has been EC-type examined by the following notified body: **CERTOTTICA S.C.R.L.**, Z.I. Villanova 32013 Longarone (BL) Italy, (Code 0530)

To download the declaration of conformity please visit: www.thigmbh.at

MARKING DEFINITIONS - THI EN166 3 F

CODE	DEFINITION
THI	Code of Manufacturer - THI GmbH
EN166	European Standard - Personal Eye Protection Specification - EN 166:2001
3	Field of use - against splashes of liquids
F	Mechanical Strength - resistance against higspeed particles - low energy 45 m/s

INDICATIONS FOR USE

The **ViVi® Helmets** are surgical devices that pull in and circulate fresh air throughout the helmet to aid user comfort. **ViVi® Helmets** are multiple-use devices, which are worn by healthcare professionals in the operating room during surgery.

CONTRAINDICATIONS

None known.

LIMITATIONS OF USE

WARNINGS:

- The surgical helmets are not intended to be used as a respiratory protection device.
- The materials contained within this device are not known allergens. However, if an allergic reaction occurs, contact a physician.

FOR USE WITH

- Rechargeable Battery Pack and Battery Holster. See *Accessories* section

DESCRIPTION

In addition to pulling in, pushing out and circulating fresh air throughout the surgical helmet, the helmet stabilizes the face-shield of a surgical gown for user comfort.



USER/PATIENT SAFETY

WARNINGS:

- Before using any component, or any component compatible with this equipment, read and understand the instructions. Pay special attention to **WARNING** information. Become familiar with the components prior to use.
- Only trained and experienced healthcare professionals should use this equipment.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage.
- DO NOT use any equipment if damage is apparent or the inspection criteria are not met. See the Inspection and Testing section for inspection criteria.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. THI GmbH, as a manufacturer, does not recommend surgical procedure or technique.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like the **ViVi® Helmet**. Install and place the helmet into service according to the EMC information in this manual. Portable and mobile RF communications equipment, such as wireless phones, can affect the function of the helmet. See the Specifications section.


SYMBOL DEFINITIONS

The symbols located on the equipment and/or labeling are defined in this section.


SYMBOL	DEFINITION
	Caution, consult accompanying documents
	Waste Electrical and Electronic Equipment Directive (WEEE)
	Recyclable Product
	Headlight ON/OFF
	Temperature Limitation
	Humidity Limitation
	Consult Instructions for Use
	Direct Current (DC)
	Hook-and-loop attachment points
	Acoustic Signal
	Manufacturer
	Non-ionizing electromagnetic radiation
	CE Mark
	Catalog number
	Serial number
	Batch number
	Quantity
	Non Sterile
	Use only in enclosed spaces (indoors)
	China RoHS
	Compliance to California Energy Commission requirements

FEATURES - HELMET


FRONT REF 80100/80120

FRONT REF 80130

BACK REF 80100, 80120 AND 80130

1	PRIMARY FAN INTAKE - The primary fan pulls in air and circulates the air throughout the helmet.
2	WASTE AIR FAN INTAKE - The secondary fan pulls in waste air from the helmet and pushes it out the back.
3	AIR DUCT - This vent channels fresh air to the front of the helmet.
4	WASTE AIR NOZZLE - This vent pushes waste air through the filter out the back of the helmet.
5	HOOK - This hook aids in positioning the hood or toga lens on the helmet.
6	ATTACHMENT POINTS (3 hook-and-loop fastener) - These fasteners secure the disposable hood or toga lens to the helmet.
7	HEIGHT ADJUSTMENT STRAP - this strap is adjusted to provide a comfortable fit. Push the release button and pull/push the strap to extend/retract the length of the strap for fit and comfort.
8	COMFORT/SOFT PADS - These removable pads provide a comfort fit.
9	HEADBAND ADJUSTMENT KNOB - Customize the fit by turning the knob clockwise to tighten and counter clockwise to loosen the headband.
10	HEADLIGHT - Available on the headlight helmet model only.
11	LIGHT SHROUD - Available on the headlight helmet model only. The shroud prevents glare from the shield.
12	HEADLIGHT ADJUSTMENT HANDLE - Available on the headlight helmet model only. The handle is used to direct the headlight beam.
13	HEADLIGHT CONTROL BUTTON - Available on the headlight helmet model only. Push the button to switch the headlight on/off.
14	POWER CABLE - This built in power cable is used to connect all helmet models with the battery pack. See the Accessories section.

INSPECTION AND TESTING

WARNINGS:

- Only individuals trained and experienced in the maintenance of reusable medical devices should inspect and test this equipment.
- DO NOT use any equipment if damage is apparent or the inspection criteria are not met.
- DO NOT disassemble or service this equipment, unless otherwise specified.

NOTE:

- For service, contact your distributor.

INTERVAL	INSPECTION CRITERIA	ACTION
Before Use	Inspect equipment for damage, wear, or missing components.	If damage is apparent, replace the equipment.
	Inspect power cord for cuts or punctures.	
	Inspect headlight lens for scratches.	
	Inspect the battery pack housing for cracks.	
	Check fan motors for excessive noise.	

NOTE: If any component must be discarded, see the *Disposal/Recycle* section.

STORAGE AND HANDLING

CAUTION:

- ALWAYS store the equipment within the specified environmental condition values throughout its useful life. See the *Specifications* section.
- ALWAYS store the battery pack between operations in the battery charger.
- To ensure the longevity, performance and safety of this equipment, use of the original packaging material is recommended when storing or transporting this equipment.

INSTRUCTIONS

WARNINGS:

- Perform recommended inspection and testing as indicated in the Inspection and Testing section.
- ALWAYS operate the equipment within the specified environmental condition values. See the Specifications section.
- ALWAYS wear the surgical helmet with a surgical toga or surgical hood to obtain the required personal protection.

CAUTIONS:

- DO NOT throw, kick, bend or twist the helmet or its features. The lightweight design is not designed to take excessive loads without damage.
- DO NOT suspend the helmet from a coat hook by using the frame of the face shield as loop.
- Do not scratch the headlight lens.
- Do not cut, stretch or kink the power cable. Failure to comply will damage the cable and may cause malfunction.

NOTE:

- Helmet position will vary from person to person depending on head size and personal comfort preferences.

TO DON A SURGICAL HELMET

WARNINGS:

- DO NOT connect the power cable plug to the battery pack before donning the helmet.

NOTE:

- ALWAYS perform the following steps in sequence.

1. Apply the comfort- or soft pads inside the helmet as required (fig. 1).

NOTE:

- Comfort- or soft pads are designed with recessions along the skin-contacting areas to promote cooling and moisture transfer. Apply the comfort- or soft pads with recessions facing towards your skin.
2. Loosen the headband with the headband adjustment knob as required by turning the knob counter clockwise (fig. 2).
 3. Adjust the helmet height using the height adjustment strap. Push the release button and pull/push the strap to extend/retract the length of the strap for fit and comfort (fig. 3).

CAUTION:

- DO NOT pull to extend or push to retract the height adjustment strap without pushing the release button. Failure to comply may result in product damage.
4. Place the helmet on your head and adjust the headband and height strap for comfort (fig. 4).

NOTES:

- Make sure the helmet fits securely on top of your head.
 - Make sure the face mask ties are accessible to allow easy removal of the face mask when donning the hood or toga.
5. Plug the power cord into the battery pack or the battery holster, once connected the fans will automatically turn on (fig. 5).

WARNINGS:

- Keep your fingers away from the fans during operation.
6. Secure the power cord into the battery pack/battery holster by tightening the connection at the base of the power cord (fig. 6).
 7. Connect the battery pack or battery holster to your waistband or put the battery into your pocket (fig. 7).
 8. Adjust the fan speed for optimal air flow by tilting your head straight back then back to your starting position. When done correctly, the helmet will emit a sound signal and the fan speed will change (fig. 8).

NOTES:

- Familiarize yourself with the gesture controlled fan speed adjustment before first time use in the OR.
- You are now ready to scrub for surgery and don a hood or toga. See the instructions for use supplied with the hood or toga.
- If the face shield lens begins to fog during use, see the Troubleshooting section.

TO USE A SURGICAL HELMET WITH A HEADLIGHT

9. Push the headlight button to switch the headlight on/off (fig. 9).

WARNINGS:

- Keep your fingers away from the headlight body during operation.
 - Avoid looking straight into the light beam.
10. Adjust the headlight beam direction with the headlight adjustment handle (fig. 10).

TO CONTROL FAN SPEED

All **Vivi®** helmets provide the user with two options to control the fan speed, manual and automatic. To manually change the fan speed the user will perform the following:

1. From a natural starting position the user should tilt their head straight back.
2. Once a "beep" sound is heard, the fan speed has changed.
3. The user returns their head to the starting position.

The "beep" sound will correspond with the speed of the fan. A low pitch sound will be heard for the low speed, two "beeps" will be heard for the middle speed and a higher pitch sound will be heard for the fastest fan speed.

The second option for the fan speed control is an automatic system built into all **Vivi®** helmets. When the helmet is powered on it starts in automatic mode; it can also be placed in automatic mode using the head gestures to reach the middle speed (two "beeps"). The automatic mode of the helmet will adjust the speed of the fan based on the activity of the user. As the user becomes more physically active, the fan will automatically begin to speed up; as activity levels begin to lower the fan speed will also lower.

NOTES:

- The user can manually switch out of automatic mode at any time using the previously described head motion.
- Familiarization with the fan speed control system is recommended to ensure user comfort during operation.

TO DOFF A SURGICAL HELMET

CAUTION:

- DO NOT doff the helmet together with the hood/toga and pull on the hood to remove it from the helmet. This might damage the helmet.
1. Loosen the power cord from the battery by turning the base of the power cord in a counter clockwise direction.
 2. Unplug the power cord from the battery, the fans will automatically turn off.
 3. Loosen the headband adjustment knob and remove the helmet from your head.
 4. Remove and discard the comfort/soft pads from the helmet. Clean the helmet and cord. See the Cleaning section. To clean the battery pack, see the instructions for use supplied with the THI battery charger.

CLEANING

WARNINGS:

- ALWAYS remove the battery pack, comfort- or soft pads from the helmet before cleaning.
- DO NOT immerse any component in liquid.
- DO NOT sterilize any component.

CAUTIONS:

- DO NOT use solvents, lubricants, or other chemicals, unless otherwise specified.
- DO NOT allow liquids or moisture into the fan or electrical connections of the helmet.
- DO NOT allow debris to enter the helmet fans.

NOTES:

- To clean the battery pack or the battery-holster, see the instructions for use supplied with the battery charger.

TO CLEAN THE SURGICAL HELMET

1. Remove all the accessories from the helmet. Discard the comfort- or soft pads.
2. Wipe the helmet with a soft cloth dampened with a nonabrasive, hospital disinfectant.
3. Dry the equipment with a lint-free towel or medical-grade compressed air.
4. Clean the headlight lens with lens tissue.
5. See the *Storage and Handling* section.

DISPOSAL/RECYCLING

WARNINGS:

• ALWAYS follow the current local recommendations and/or regulations governing environmental protection and risks associated with recycling or disposing of the equipment at the end of its useful life.

- ALWAYS decontaminate equipment exposed to infectious material before sending the equipment to a waste treatment facility.
- Per the European Union (EU) Waste Electronics and Electrical Equipment (WEEE) Directive 2012/19/EU, product must be collected separately.
- DO NOT dispose of as unsorted municipal waste. Contact your local distributor for disposal information.
- Contains no hazardous substances, per RoHS standard.

ACCESSORIES

This section describes components that may be ordered to replace original equipment that is damaged, worn, or must be replaced. This section may also contain optional components used with the surgical helmet.

WARNINGS:

- Use only THI GmbH-approved components and accessories.
- Using other electronic components and accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the equipment.
- DO NOT modify any component or accessory.
- DO NOT reuse, reprocess, or repackage a single use device. A single use device is intended for a single use only. A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the single use device resulting in operational failure. Critical product information may be lost if the single use device is re-packaged. Failure to comply may lead to infection or cross-infection and result in patient and/or health care staff injury.
- Make sure the power plug is easily accessible when connected to the mains.

The following THI GmbH-approved accessories are sold separately:

DESCRIPTION	REF
Rechargeable Battery Pack	60401
Rechargeable Battery Pack HW	60402
Battery Holster HW	60403
Vivi® Comfort Pads – A set of disposable pads are installed in the inner helmet to provide additional comfort; one set per case.	50785
Vivi® Soft Pads – A set of disposable pads are installed in the inner helmet to provide extra soft comfort; one set per case.	50786
NOTE: For a complete list of accessories, contact your distributor.	

TROUBLESHOOTING


WARNING:

- DO NOT service this equipment, unless otherwise specified. If you require service, contact your local distributor.

PROBLEM	CAUSE	ACTION
The fans fail to provide air flow.	The battery pack is not charged or is damaged.	Change or replace the battery pack.
	The battery is not correctly connected.	Check the battery connection.
	The fans are defective or the connection is broken.	Replace the helmet.
The fans fail to turn on when connected to the battery pack.	The battery requires replacement.	Replace the battery pack with a fully charged battery pack.
	The power cord is damaged or the fans are damaged.	Replace the helmet.
	The battery is not correctly connected.	Check the battery connection.
The face shield lens is fogging.	Not enough airflow is circulating throughout the helmet.	Use the required head motion to increase speed of the fans.
	The battery pack requires replacement.	Replace the battery pack with a fully charged battery pack.
	The helmet is damaged.	Replace the helmet.
The helmet feels unbalanced.	The headbelt is not properly adjusted.	Tighten or loosen the headbelt using the headbelt adjustment knob.
	The comfort/soft pads are not properly placed.	Correct placement of the comfort pads. See <i>Instructions</i> section.
	The height strap is not properly adjusted.	Raise or lower height strap on the helmet. See <i>Instructions</i> section.
	The helmet is damaged.	Replace the helmet.
The head gesture does not change the fan speed.	The helmet is not correctly placed on the users head.	Adjust the placement of the helmet on the users head.
	The helmet is damaged	Replace the helmet
	The helmet is not correctly placed on the users head.	Adjust the placement of the helmet on the users head. See <i>Fan Speed Control</i> section.
Sporadic electrical interference is experienced.	The head gesture is done incorrectly.	See <i>Fan Speed Control</i> section.
	Electrical noise is present.	Turn off all the electrical equipment not in use in the room, including cell phones. Relocate the electrical equipment to maximize the distance between the equipment. Increase spatial distance. Plug equipment into different outlets.

NOTE: If any component must be discarded, see the *Disposal/Recycle* section

TECHNICAL SPECIFICATIONS - HELMETS

MODEL	REF 80100	REF 80120	REF 80130
Description	viVi® Helmet	viVi® Helmet with Headlight	viVi® Helmet HPL
European Conformity	CE		
Mass	0.37 kg [13.1 oz]	0.41 kg [14.5 oz]	0.40 kg [14.2 oz]
Fan Speed	Primary Fan: 3,925 - 4,775 rpm, Waste Air Fan: 5,885 rpm		
Input Voltage	12.4 - 16.4 V  (DC)		
Ingress Protection (IP)	IPX0		
Mode of Operation	Continuous		

ENVIRONMENTAL CONDITIONS	OPERATION	STORAGE AND TRANSPORTATION
Temperature Limitation	min. 10° C  max. 40° C	min. -20° C  max. 40° C
Humidity Limitationf	min. 30 %  max. 75 %	min. 10 %  max. 75 %

SPECIFICATIONS

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

viVi® Helmet REF 80100, 80120 and 80130 are intended for use in the electromagnetic environment specified below.
The customer or the user of viVi® Helmet REF 80100, 80120 and 80130 should assure that it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	viVi® Helmet REF 80100, 80120 and 80130 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class A	
Harmonic emissions IEG 61000-3-2	Not Applicable	viVi® Helmet REF 80100, 80120 and 80130 are suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEG 61000-3-3	Not Applicable	

CHARGING INSTRUCTIONS

WARNING:

ALWAYS use the REF 60500 Battery Charger to charge and store REF 60401 battery pack and REF 60501 single Battery Charger HW for REF 60402 battery pack.

CAUTION:

For tips and recommendations on how to handle the Battery Pack REF 60401/ REF 60402 refer to the enclosed operating instructions and the dedicated operation instructions of battery pack and battery charger.

NOTE:

For optimum performance, charge the battery before first use and after each use. Check charging status by inserting the battery into the charging bay of the battery charger and observing the indication lights. Battery Pack REF 60402 has charging information display (Pos. 5) included.

- Place the battery charger on a flat surface away from sources of heat and humidity.
- Connect the DC connector on the power supply (Pos. 4) with the battery charger.

CAUTION:

DO NOT connect power adaptors other than provided by THI to your battery charger.

- Plug the power supply (REF 60500), using the AC head for your power system or the power cord (REF 60501) into the mains outlet.
- a) - for REF 60500: Check if the power LED (Pos. 2) on the power adaptor and the power LED (Pos. 11) on the battery charger are ON.
b) - for REF 60501: Right after plugging in the power cord (Pos. 1) the charging LED (Pos. 10) is showing red-yellow-green, then the LED is going out (charger self test)
- Insert a battery pack (Pos. 7) into the charging bay (Pos. 6)

NOTE:

Make sure that the battery is orientated correctly in the charging bay. Only if the charging contacts (Pos. 9) are in-line with its counterparts in the charging bay, a connection is established.

- The charging indication (Pos. 10) will indicate the charging status (see charging indication) and the charger will automatically charge the battery.
- After charging is completed, the charger will maintain the battery pack in a fully charged state.

NOTE:

We recommend leaving the battery pack in the charging bay if they are intended for use.

TROUBLESHOOTING

WARNING:

DO NOT service this equipment, unless otherwise specified. If you require service, contact your local distributor.

PROBLEM	CAUSE	ACTION
Battery not charging	No power	Check power indications on charger and power adaptor. Make sure the power adaptor is connected to the AC outlet and the charger.
	Incorrect battery placement	Make sure the battery is well-seated in the charging bay.
	Battery defective	Replace battery.
	Battery charger or power adaptor defective	Replace equipment.
'Error' indication at individual charging bay (red LED)	Battery 'overcharged'	Discharge/charge the battery once.
	Defective battery	Replace battery.

NOTE: If any component must be discarded, see the *Disposal/Recycle* section

BATTERY DISPOSAL/RECYCLING

WARNING:

• ALWAYS follow the current local recommendations and/or regulations governing environmental protection and risks associated with recycling or disposing of the equipment at the end of its useful life.

- Per the European Union (EU) Waste Electronics and Electrical Equipment (WEEE) Directive 2012/19/EU, product must be collected separately.
- DO NOT dispose of as unsorted municipal waste. Contact your local distributor for disposal information.
- Contains no hazardous substances, per RoHS standard.

NOTES:

Li-ion

- This battery pack contains lithium ionide. After the end of its useful life, the battery must be disposed of by qualified recycling or hazardous materials handling operators.
- DO NOT dispose the battery pack between other solid wastes. For disposal and recycling information, refer to local operators responsible for recycling or handling hazardous materials.

TECHNICAL SPECIFICATION - BATTERY

MODEL	REF 60401	REF 60402
Description:	Rechargeable Battery Pack	
European Conformity:	CE	
Cell chemistry:	Lithium Ion	
Voltage (nominal):	14.8 V (DC)	14.4 V (DC)
Capacity (rated):	2.950 mAh (43.6 Wh)	3.450 mAh (49.7 Wh)
Cycle life characteristics:	> 80% of rated capacity after 400 duty cycles	> 75% of rated capacity after 300 duty cycles
Ingress Protection (IP):	IPX0	
Mode of operation:	Continuous	
Charging time:	4 hours	2 ½ hours
Battery life:	> 6 hours (with REF 80100 vivi® Helmet) > 3 ½ hours (with REF 80130 vivi® Helmet HPL)	> 7 hours (with REF 80100 vivi® Helmet) > 4 hours (with REF 80130 vivi® Helmet HPL)
Weight:	250 g (8.8 oz)	230 g (8.1 oz)


CHARGING INDICATION

NOTE:

The battery status is indicated by a multi-colour LED (Pos.11) at charging bay.


LED	STATUS INFORMATION REF 60500	STATUS INFORMATION REF 60501
OFF	No battery in bay or no power.	No battery in bay or no power.
Solid green	Stand by/battery present/battery full	Stand by/battery present/battery full
Solid orange	Battery charging	Battery charging
Solid red	Error (see troubleshooting section)	Either: • The battery is badly damaged and needs to be replaced or • It is a conventional battery which can not be charged
Red blinking	- - -	The battery is too hot or too cold to be charged without damage

FEATURES - BATTERY AND CHARGER REF 60500

1	AC/DC POWER ADAPTOR - with multiple AC Heads.
2	POWER LED - A blue LED on the power adaptor indicating power ON.
3	AC HEADS - 4 different power plugs for country specific use. Consult the user guide of the power adaptor how to install a suitable AC head.
4	ADAPTOR CORD - Connect the DC connector on the adaptor cord with the battery charger.
5	BATTERY CHARGER - designed to charge 4 battery packs REF 60401 simultaneously. Use the charger to store battery packs between operations.
6	CHARGING BAY - Insert the battery pack into the charging bay to connect the battery with the charger. Ensure correct orientation of the battery.
7	BATTERY PACK - a belt clip helps to secure the battery pack on the body.
8	BATTERY CONNECTOR - This built in connector is used to connect the battery with the power cord of all VIVI® helmets. Consult the section To Don a Surgical Helmet how to secure the power cord into the battery pack.
9	CHARGING CONTACTS - Ensure they are dry and clean when connecting.
10	CHARGING INDICATION - Each charging bay is equipped with a multi-color LED light indicating charging status. Fully charged Charging Error
11	POWER LED - A green LED on the battery charger indicating power ON.

FEATURES - BATTERY AND CHARGER REF 60501

1	AC/DC POWER SUPPLY with power cord.
2	Battery Holster HW - The battery holster is used to connect the battery with the power cable of the helmet. Insert the battery pack in the battery holster. Ensure correct orientation of the battery.
3	BATTERY CHARGER HW - designed to one battery packs REF 60402. Use the charger to store battery packs between operations.
4	POWER SUPPLY CORD - Connect the DC connector on the Power Supply cord with the battery charger.
5	CHARGING INFORMATION DISPLAY - Displays information about the charging state of the battery by pressing the circle symbol.
6	CHARGING BAY - Insert the battery pack into the charging bay to connect the battery with the charger. Ensure correct orientation of the battery.
7	BATTERY PACK
8	BATTERY CONNECTOR - This built in connector is used to connect the battery with the power cord of all VIVI® helmets. Consult the section To Don a Surgical Helmet how to secure the power cord into the battery pack.
9	CHARGING CONTACTS - Ensure they are dry and clean when connecting.
10	CHARGING INDICATION - The charging bay is equipped with a multi-color LED light indicating charging status.

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